

Safety Blood Collection Set

Secure-Lock Type (with Luer Adapter)

Product Advantage:

- ISO color coding for easy recognition
- Luer Lock and Luer slip adapter available
- Quality needles for patient comfort
- One-hand operation
- Award-winning secure-locking safety design to prevent needle-stick injury and reuse



Safety Secure:

- After blood collecting, press the the cap hard until a click is heard;
- The needle will be locked in the sheath permanently;

The risks from needle-stick injuries can largely decrease with other unsafe injection practices, such as needle reuse and improper medical waste collection.

Product Specification:

Gauge	OD(mm)	Color	Length
23G	0.60	blue ●	1" ~ 1 1/2" (25~38mm)
22G	0.70	black ●	1" ~ 1 1/2" (25~38mm)
21G	0.80	green ●	1" ~ 1 1/2" (25~38mm)
20G	0.90	yellow ●	1" ~ 1 1/2" (25~38mm)
18G	1.20	pink ●	1" ~ 1 1/2" (25~38mm)

