

Safety Blood Collection Set -Pen Type (with or without tube holder)

Product Advantage:

- Safe:
 - * Safety cap prevent needle-stick injuries and needle reuse;
 - * Top protect cap and bottom protect cap double secure the hygiene and sterilization.
- Convenient:
 - * Activate the safety cap only with one hand.
- Sharp bevel, less pain.



Safety Secure:

- After blood collecting, press the the cap harder until a click is heard;
- The needle will be locked in the sheath permanently;

The risks from needle-stick injuries can largely decrease with other unsafe injection practices, such as needle reuse and improper medical waste collection.



Product Specification:

* Below is the normal size, we can also provide other size upon client's request.

Gauge	OD(mm)	Color	Length
23G	0.60	blue ●	1"~1 1/2" (25~38mm)
22G	0.70	black ●	1"~1 1/2" (25~38mm)
21G	0.80	green ●	1"~1 1/2" (25~38mm)
20G	0.90	yellow ●	1"~1 1/2" (25~38mm)
18G	1.20	pink ●	1"~1 1/2" (25~38mm)